

Anti-Ire1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O75460

Item Number: CM0024684



Introduction

Rabbit polyclonal antibody against human IRE1 (ERN1). Validated for WB, IF/ICC, IHC-P, ELISA. Detects IRE1 in human, mouse, rat. Key sensor of ER stress and UPR pathway.

[Learn More](#)

Parameter	Specification
Product Name	IRE1 Rabbit pAb
Remarks/Alias	IRE1; IRE1P; IRE1a; hIRE1p
Species	Human
GeneID (Human)	2081
GeneID	2081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 286-435 of human IRE1 (NP_001424.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75460
Protein Weight	110 kDa
Shipping	Ice bag