

Anti-Irf1 Rabbit Polyclonal Antibody For Wb, Elisa - P15314

Item Number: CM0028540



Introduction

Rabbit polyclonal antibody targeting IRF1, validated for WB and ELISA with cross-reactivity to mouse. Suitable for studying interferon regulatory factor 1 functions in immunity, inflammation, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	IRF1 Rabbit pAb
Remarks/Alias	Irf-1
Species	Mouse
Gene ID (Human)	3659
Gene ID	16362
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 113-212 of human IRF1 (NP_032416.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P15314
Protein Weight	37 kDa
Shipping	Ice bag