


Anti-Irf2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - P14316

Item Number: CM0011953



Introduction

High-quality IRF2 Rabbit Monoclonal Antibody validated for WB, IP, ELISA. Specific for human and mouse IRF2 (P14316). Ideal for transcription regulation and interferon signaling research.

[Learn More](#)

Parameter	Value
Product Name	IRF2 Rabbit mAb
Remarks/Alias	IRF-2; IRF2
Species	Human
GeneID (Human)	3660
GeneID	3660
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P14316
Protein Weight	39kDa
Shipping	Ice bag