


Anti-Irf2Bp2 Polyclonal Antibody For Wb, Elisa - Q7Z5L9

Item Number: CM0020887



Introduction

High-quality rabbit polyclonal antibody targeting human IRF2BP2 (Q7Z5L9). Validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for transcriptional regulation and B-cell differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	IRF2BP2 Rabbit pAb
Remarks/Alias	CVID14; IRF2BP2
Species	Human
GeneID (Human)	359948
GeneID	359948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 475-545 of human IRF2BP2 (Q7Z5L9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z5L9
Protein Weight	61kDa
Shipping	Ice bag