

Anti-Irf3 (Iiae7) Rabbit Polyclonal Antibody For Wb And Elisa - Q764M6

Item Number: CM0028542

Introduction

IRF3 Rabbit Polyclonal Antibody (IIAE7), validated for WB and ELISA detection, targets Sus scrofa IRF3 with human cross-reactivity. Suited for innate immune and antiviral studies.

[Learn More](#)

Parameter	Value
Product Name	IRF3 Rabbit pAb
Remarks/Alias	IIAE7
Species	Sus scrofa
GeneID (human)	3661
GeneID	396656
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of Sus scrofa IRF3 (NP_998935.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q764M6
Protein Weight	47kDa
Shipping	Ice bag