

Anti-Irf3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - Q14653

Item Number: CM0013362



Introduction

IRF3 Rabbit Polyclonal Antibody (pAb) for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Validated in human, mouse, rat. Targets interferon regulatory factor 3 (Q14653), key innate immune transcription factor. Ideal for antiviral and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	IRF3 Rabbit pAb
Remarks/Alias	IIAE7; IRF3
Species	Human
Gene ID (Human)	3661
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14653
Protein Weight	47kDa
Shipping	Ice bag