

Anti-Irf4 Polyclonal Antibody For Wb, Elisa - Q15306

Item Number: CM0017229



Introduction

Rabbit polyclonal anti-IRF4 antibody validated for Western blotting and ELISA, recognizing human, mouse, and rat IRF4. Suitable for studying IRF4 transcriptional activation in lymphoid and immune cells.

[Learn More](#)

Parameter	Value
Product Name	IRF4 Rabbit pAb
Remarks/Alias	MUM1; LSIRF; SHEP8; NF-EM5; IRF4
Species	Human
Gene ID (Human)	3662
Gene ID	3662
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-350 of human IRF4 (NP_002451.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15306
Protein Weight	52kDa
Shipping	Ice bag