


Anti-Irf5 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q13568

Item Number: CM0002376



Introduction

Rabbit monoclonal antibody targeting human IRF5 (SLEB10), validated for WB, IP, ELISA. Cross-reacts with mouse. Ideal for studying innate immunity and interferon signaling.

[Learn More](#)

Parameter	Value
Product Name	IRF5 Rabbit mAb
Remarks/Alias	SLEB10; IRF5
Species	Human
Gene ID (Human)	3663
Gene ID	3663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 6-151 of human IRF5 (Q13568)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13568
Protein Weight	56kDa
Shipping	Ice bag