


Anti-Irf5 Rabbit Polyclonal Antibody For Wb - Q13568

Item Number: CM0015230



Introduction

Rabbit polyclonal antibody against human IRF5, validated for WB, IHC-P, and ELISA. Suitable for studies on innate immunity, interferon signaling, and inflammatory response. Cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	IRF5 Rabbit pAb
Remarks/Alias	SLEB10; IRF5
Species	Human
GeneID (human)	3663
GeneID	3663
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human IRF5 (NP_001092097.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13568
Protein Weight	56kDa
Shipping	Ice bag