

Anti-Irf6 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - O14896

Item Number: CM0010228

Introduction

IRF6 Rabbit Monoclonal Antibody (O14896) validated for WB, IF/ICC, IP, ELISA. Detects human, mouse, rat IRF6. Ideal for epidermal differentiation and breast cancer research.

[Learn More](#)

Parameter	Specification
Product Name	IRF6 Rabbit mAb
Remarks/Alias	LPS; PIT; PPS; VWS; OFC6; PPS1; VWS1; IRF6
Species	Human
GeneID (Human)	3664
GeneID	3664
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14896
Protein Weight	53kDa
Shipping	Ice bag