

Anti-Irf7 Rabbit Polyclonal Antibody For Wb, Elisa - Q92985

Item Number: CM0014566



Introduction

Rabbit polyclonal antibody against human IRF7 for WB and ELISA; cross-reacts with mouse. Key transcriptional regulator of type I interferon response, critical for innate immunity against viruses.

[Learn More](#)

Parameter	Specification
Product Name	IRF7 Rabbit pAb
Remarks/Alias	IMD39; IRF-7; IRF7A; IRF7B; IRF7C; IRF7H; IRF-7H; IRF7
Species	Human
GeneID (Human)	3665
GeneID	3665
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human IRF7 (NP_001563.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q92985
Protein Weight	54kDa
Shipping	Ice bag