

Anti-Irf8 Polyclonal Antibody For Wb, Elisa - Q02556

Item Number: CM0015966

Introduction

Rabbit polyclonal antibody against human IRF8 (ICSBP), validated for Western Blot and ELISA. IRF8 is a transcription factor involved in interferon signaling and dendritic cell development. Shipped on ice.

[Learn More](#)

Parameter	Value
Product Name	IRF8 Rabbit pAb
Remarks/Alias	ICSBP; IRF-8; ICSBP1; IMD32A; IMD32B; H-ICSBP; IRF8
Species	Human
Gene ID (Human)	3394
Gene ID	3394
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 137-426 of human IRF8 (NP_002154.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q02556
Protein Weight	48kDa
Shipping	Ice bag