

Anti-Irf8 Rabbit Polyclonal Antibody For Wb And Elisa - Q02556

Item Number: CM0028355

Introduction

Rabbit polyclonal antibody against IRF8, validated for WB and ELISA. Reacts with human and mouse. Ideal for studying transcription factor involved in immune cell differentiation and interferon signaling.

[Learn More](#)

Parameter	Value
Product Name	IRF8 Rabbit pAb
Remarks/Alias	ICSBP, IRF-8, ICSBP1, IMD32A, IMD32B, H-ICSBP, IRF8
Species	Human
GeneID (human)	3394
GeneID	3394
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 137-426 of human IRF8 (NP_002154.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q02556
Protein Weight	48kDa
Shipping	Ice bag