

Anti-Irp2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P48200

Item Number: CM0004234



Introduction

High-quality Anti-IRP2 Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Shows cross-reactivity with human, mouse, and rat. Ideal for iron metabolism, gene regulation, and disease biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	IRP2 Rabbit mAb
Remarks/Alias	ACO3; IRP2; IRP2AD; NDCAMA; IRE-BP2; IRE-BP 2
Species	Human
Gene ID (Human)	3658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-963 of human IRP2 (NP_004127.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48200
Protein Weight	105 kDa
Shipping	Ice bag