

Anti-Irp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P48200

Item Number: CM0028818



Introduction

Rabbit polyclonal antibody against human IRP2 (IREB2), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for iron metabolism and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	IRP2 Rabbit pAb
Remarks/Alias	ACO3; IRP2; IRP2AD; NDCAMA; IRE-BP2; IRE-BP 2
Species	Human
Gene ID (Human)	3658
Gene ID	3658
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 760-963 of human IRP2 (NP_004127.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48200
Protein Weight	105kDa
Shipping	Ice bag