

Anti-Irs1 Rabbit Polyclonal Antibody For Wb And Elisa - P35569

Item Number: CM0027926



Introduction

Rabbit polyclonal antibody against mouse IRS1 (insulin receptor substrate 1), suitable for Western blot and ELISA. Useful for research on insulin signaling, glucose metabolism, and cell growth pathways.

[Learn More](#)

Parameter	Value
Product Name	IRS1 Rabbit pAb
Remarks/Alias	G972R; IRS-1
Species	Mouse
Gene ID (Human)	3667
Gene ID	16367
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P35569
Protein Weight	131 kDa
Shipping	Ice bag