

Anti-Irs1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P35568

Item Number: CM0014028



Introduction

Rabbit polyclonal antibody against human IRS1, validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat IRS1. Ideal for insulin signaling and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	IRS1 Rabbit pAb
Remarks/Alias	HIRS-1; IRS1
Species	Human
Gene ID (Human)	3667
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1043-1242 of human IRS1 (NP_005535.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35568
Protein Weight	132kDa
Shipping	Ice bag