

Anti-Irs2 Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Y4H2

Item Number: CM0008823



Introduction

Rabbit monoclonal antibody against human IRS2, suitable for Western Blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Predicted molecular weight 137 kDa. Ideal for insulin signaling and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	IRS2 Rabbit mAb
Remarks/Alias	IRS-2, Irs2, IRS 2
Species	Human
Gene ID (Human)	8660
Gene ID	8660
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1239-1338 of human IRS2 (NP_003740.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4H2
Protein Weight	137kDa
Shipping	Ice bag