

Anti-Iscu Rabbit Polyclonal Antibody For Wb, Elisa - Q9H1K1

Item Number: CM0026244



Introduction

Rabbit polyclonal antibody raised against recombinant human ISCU protein, suitable for WB and ELISA. Validated for human samples; targets the mitochondrial scaffold protein involved in iron-sulfur cluster biogenesis.

[Learn More](#)

Parameter	Value
Product Name	ISCU Rabbit pAb
Remarks/Alias	HML; ISU2; NIFU; NIFUN; hnifU; 2310020H20Rik; ISCU
Species	Human
GeneID (human)	23479
GeneID	23479
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-167 of human ISCU (NP_998760.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H1K1
Protein Weight	18kDa
Shipping	Ice bag