

Anti-Isg20 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Az6

Item Number: CM0015524



Introduction

Rabbit polyclonal antibody against human ISG20, validated for WB and ELISA. Cross-reacts with mouse and rat. Useful for studies on antiviral innate immunity, interferon-stimulated genes, and RNA metabolism.

[Learn More](#)

Parameter	Specification
Product Name	ISG20 Rabbit pAb
Remarks/Alias	CD25, HEM45, ISG20
Species	Human
Gene ID (Human)	3669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human ISG20 (NP_002192.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AZ6
Protein Weight	20kDa
Shipping	Ice bag