

Anti-Islet1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P61371

Item Number: CM0002219

Introduction

Rabbit monoclonal antibody targeting human Islet1 (P61371), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. Ideal for studying neuronal differentiation and heart development.

[Learn More](#)

Parameter	Value
Product Name	Islet1 Rabbit mAb
Remarks/Alias	Isl-1; ISLET1; Islet1
Species	Human
GeneID (Human)	3670
GeneID	3670
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Islet1 (P61371).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P61371
Protein Weight	39 kDa
Shipping	Ice bag