

Anti-Islet1 Rabbit Polyclonal Antibody For Wb, Elisa - P61371

Item Number: CM0028848

Introduction

High-quality anti-Islet1 rabbit polyclonal antibody validated for WB and ELISA applications. Reacts with human, mouse, and rat Islet1, a key transcription factor in development. Ideal for neuroscience, cardiac, and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Islet1 Rabbit pAb
Remarks/Alias	Isl-1; ISLET1; Islet1
Species	Human
GeneID (human)	3670
GeneID	3670
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-349 of human Islet1 (NP_002193.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61371
Protein Weight	39kDa
Shipping	Ice bag