


Anti-Ist1 Polyclonal Antibody For Wb, Elisa - P53990

Item Number: CM0029754



Introduction

Rabbit polyclonal antibody against IST1 (CHMP8, OLC1). Validated for WB and ELISA, reacts with human, mouse, and rat. Suitable for studying ESCRT-III-mediated cytokinesis and nuclear envelope dynamics.

[Learn More](#)

Parameter	Value
Product Name	IST1 Rabbit pAb
Remarks/Alias	OLC1; CHMP8; IST1
Species	Human
GeneID (Human)	9798
GeneID	9798
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human IST1 (NP_001257906.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53990
Protein Weight	40kDa
Shipping	Ice bag