


Anti-Isy1 Polyclonal Antibody For Wb, Elisa - Q9ULr0

Item Number: CM0018163



Introduction

Rabbit polyclonal antibody targeting human ISY1, a pre-mRNA splicing factor and miRNA processing regulator. Validated for Western blot and ELISA applications, with cross-reactivity to human and mouse samples. Ideal for splicing and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	ISY1 Rabbit pAb
Remarks/Alias	FSAP33; ISY1
Species	Human
Gene ID (Human)	57461
Gene ID	57461
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-307 of human ISY1 (NP_001186398.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULR0
Protein Weight	33kDa
Shipping	Ice bag