

Anti-Isyna1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nph2

Item Number: CM0027086

Introduction

Rabbit polyclonal antibody targeting human ISYNA1, validated for Western blot (WB) and ELISA. Detects human, mouse, and rat ISYNA1. ISYNA1 is the rate-limiting enzyme in myo-inositol biosynthesis. Suitable for metabolic and lipid signaling research.

[Learn More](#)

Parameter	Value
Product Name	ISYNA1 Rabbit pAb
Remarks/Alias	IPS; INO1; INOS; IPS 1; IPS-1; ISYNA1
Species	Human
Gene ID (Human)	51477
Gene ID	51477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human ISYNA1 (NP_057452.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPH2
Protein Weight	61kDa
Shipping	Ice bag