


Anti-Itch Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96J02

Item Number: CM0022728



Introduction

Rabbit polyclonal antibody against ITCH (AIP4), validated for WB, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying E3 ubiquitin ligase function and inflammatory signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ITCH Rabbit pAb
Remarks/Alias	AIF4; AIP4; ADMFD; NAPP1; ITCH
Species	Human
GeneID (Human)	83737
GeneID	83737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-155 of ITCH (NP_001244066.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96J02
Protein Weight	103kDa
Shipping	Ice bag