


Anti-Itch Rabbit Monoclonal Antibody (Ko Validated) For Wb, Elisa - Q96J02

Item Number: CM0010554



Introduction

Rabbit monoclonal antibody specific to ITCH (Q96J02), an E3 ubiquitin ligase. Validated for WB, ELISA; cross-reacts with human, mouse, and rat. Suitable for studying ubiquitination and immune signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ITCH Rabbit mAb
Remarks/Alias	AIF4; AIP4; ADMFD; NAPP1; CH
Species	Human
Gene ID (Human)	83737
Gene ID	83737
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ITCH (Q96J02).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96J02
Protein Weight	103kDa
Shipping	Ice bag