

Anti-Itfg1 Polyclonal Antibody For Wb, Elisa - Q8Tb96

Item Number: CM0013338

Introduction

Rabbit polyclonal antibody against human ITFG1, suitable for WB and ELISA applications. Cross-reacts with mouse and rat. Target is a T-cell immunomodulatory protein involved in graft-versus-host disease protection.

[Learn More](#)

Parameter	Value
Product Name	ITFG1 Rabbit pAb
Remarks/Alias	TIP; CDA08; LNKN-1; 2310047C21Rik; ITFG1
Species	Human
Gene ID (Human)	81533
Gene ID	81533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 313-612 of human ITFG1 (NP_110417.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TB96
Protein Weight	68kDa
Shipping	Ice bag