

Anti-Itga10 Rabbit Polyclonal Antibody For Wb And Elisa - 075578

Item Number: CM0024281

Introduction

Rabbit polyclonal antibody against human ITGA10, validated for WB and ELISA. Reacts with human and mouse. Suitable for collagen receptor studies. Protein weight: 128kDa.

[Learn More](#)

Parameter	Value
Product Name	ITGA10 Rabbit pAb
Remarks/Alias	PRO827; ITGA10
Species	Human
GeneID (Human)	8515
GeneID	8515
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 801-1040 of human ITGA10 (NP_003628.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75578
Protein Weight	128kDa
Shipping	Ice bag