


Anti-Itga1/Cd49A Monoclonal Antibody () For Fc - P56199

Item Number: CM0012666



Introduction

-conjugated rabbit anti-human ITGA1/CD49a monoclonal antibody for flow cytometry. Recognizes integrin alpha-1, a laminin and collagen receptor, enabling reliable detection of this cell adhesion molecule on human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ITGA1/CD49a mAb
Remarks/Alias	VLA1; CD49a
Species	Human
Gene ID (Human)	3672
Gene ID	3672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of human ITGA1/CD49a (NP_852478.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P56199
Protein Weight	131kDa
Shipping	Ice bag