

Anti-Itga8 Rabbit Polyclonal Antibody For Wb And Elisa - P53708

Item Number: CM0017333

Introduction

Rabbit polyclonal antibody against ITGA8, raised with recombinant protein immunogen. Validated for WB and ELISA, cross-reactive with human, mouse, and rat. Target molecular weight 117 kDa. ITGA8 is an integrin alpha subunit involved in kidney development and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	ITGA8 Rabbit pAb
Remarks/Alias	ITGA8
Species	Human
GeneID (Human)	8516
GeneID	8516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 907-1000 of human ITGA8 (NP_003629.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53708
Protein Weight	117kDa
Shipping	Ice bag