

Anti-Itga9 Polyclonal Antibody For Wb, Elisa - Q13797

Item Number: CM0025116

Introduction

Rabbit polyclonal antibody targeting Integrin alpha-9 (ITGA9), validated for WB and ELISA with mouse cross-reactivity. Ideal for cell adhesion and vascular research studies.

[Learn More](#)

Parameter	Value
Product Name	ITGA9 Rabbit pAb
Remarks/Alias	RLC; ITGA4L; ALPHA-RLC; ITGA9
Species	Human
GeneID (human)	3680
GeneID	3680
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-250 of human ITGA9 (NP_002198.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13797
Protein Weight	114kDa
Shipping	Ice bag