

Anti-Itgae (Cd103) Rabbit Polyclonal Antibody For Wb And Elisa - P38570

Item Number: CM0026372



Introduction

High-quality rabbit polyclonal antibody targeting integrin alpha-E (CD103). Validated for Western blot, ELISA with cross-reactivity to human, mouse, and rat. Ideal for investigating lymphocyte adhesion, mucosal immunity, and epithelial barrier biology.

[Learn More](#)

Parameter	Value
Product Name	ITGAE Rabbit pAb
Remarks/Alias	CD103; HUMINAE; ITGAE
Species	Human
Gene ID (Human)	3682
Gene ID	3682
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-200 of human ITGAE (NP_002199.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38570
Protein Weight	130kDa
Shipping	Ice bag