

Anti-Itgb3Bp Polyclonal Antibody For Wb And Elisa - Q13352

Item Number: CM0030324



Introduction

Rabbit polyclonal antibody targeting human ITGB3BP (CENPR, NRIF3), validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for transcription regulation and centromere/kinetochore research.

[Learn More](#)

Parameter	Value
Product Name	ITGB3BP Rabbit pAb
Remarks/Alias	CENPR; NRIF3; TAP20; CENP-R; HSU37139; ITGB3BP
Species	Human
GeneID (Human)	23421
GeneID	23421
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ITGB3BP (NP_055103.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13352
Protein Weight	20kDa
Shipping	Ice bag