


Anti-Itgb5 Rabbit Polyclonal Antibody For Wb And Elisa - P18084

Item Number: CM0024617



Introduction

ITGB5 Rabbit Polyclonal Antibody for WB and ELISA. Detects human, mouse, rat ITGB5. Recombinant immunogen, 88 kDa. For cell adhesion and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	ITGB5 Rabbit pAb
Remarks/Alias	ITGB5
Species	Human
Gene ID (Human)	3693
Gene ID	3693
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 490-720 of human ITGB5 (NP_002204.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag