

Anti-Itgb6 Rabbit Polyclonal Antibody For Wb, Elisa - P18564

Item Number: CM0019719



Introduction

Rabbit polyclonal antibody against human ITGB6 for Western blot (WB) and ELISA. Recognizes human, mouse, and rat ITGB6. Integrin beta-6 receptor involved in cell adhesion, TGF-beta activation, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	ITGB6 Rabbit pAb
Remarks/Alias	AI1H; ITGB6
Species	Human
Gene ID	3694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 732-788 of human ITGB6 (NP_000879.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18564
Protein Weight	86kDa
Shipping	Ice bag