

Anti-Itgb7 Polyclonal Antibody For Wb, Elisa - P26010

Item Number: CM0030086



Introduction

ITGB7 rabbit polyclonal antibody for Western blot and ELISA, validated in human and mouse. Targets integrin beta-7, a key lymphocyte homing receptor for gut-associated lymphoid tissue research.

[Learn More](#)

Parameter	Value
Product Name	ITGB7 Rabbit pAb
Remarks/Alias	ITGB7
Species	Human
Gene ID (Human)	3695
Gene ID	3695
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ITGB7 (NP_000880.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26010
Protein Weight	87kDa
Shipping	Ice bag