

Anti-Itgb8 Rabbit Polyclonal Antibody For Wb And Elisa - P26012

Item Number: CM0029509



Introduction

Rabbit polyclonal antibody to human ITGB8, validated for WB and ELISA. Cross-reacts with mouse and rat. Targets integrin beta-8, a fibronectin receptor involved in TGF-beta activation and vasculogenesis.

[Learn More](#)

Parameter	Value
Product Name	ITGB8 Rabbit pAb
Remarks/Alias	ITGB8
Species	Human
Gene ID (Human)	3696
Gene ID	3696
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-320 of human ITGB8 (NP_002205.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26012
Protein Weight	86 kDa
Shipping	Ice bag