

Anti-Itgbl1 Rabbit Polyclonal Antibody For Elisa And Wb - O95965

Item Number: CM0021863

Introduction

ITGBL1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat ITGBL1. Ideal for studying integrin beta-like protein 1 function in vascular biology and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	ITGBL1 Rabbit pAb
Remarks/Alias	OSCP; TIED; ITGBL1
Species	Human
GeneID (human)	9358
GeneID	9358
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human ITGBL1 (NP_004782.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95965
Protein Weight	54kDa
Shipping	Ice bag