

Anti-Itih3 Rabbit Polyclonal Antibody For Wb, Elisa - Q06033

Item Number: CM0025451

Introduction

Rabbit polyclonal antibody against human ITIH3, suitable for WB and ELISA. Cross-reacts with mouse and rat. Detects a ~100kDa band. For hyaluronan metabolism and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	ITIH3 Rabbit pAb
Remarks/Alias	H3P; SHAP; ITI-HC3; ITIH3
Species	Human
GeneID (human)	3699
GeneID	3699
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-180 of human ITIH3 (NP_002208.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06033
Protein Weight	100kDa
Shipping	Ice bag