

Anti-Itk/Emt Monoclonal Antibody For Wb, Elisa - Q08881

Item Number: CM0005451



Introduction

Rabbit monoclonal antibody targeting human ITK/EMT (tyrosine-protein kinase ITK/TSK), validated for Western blot and ELISA. Essential for T-cell signaling, immune response, and kinase activity studies.

[Learn More](#)

Parameter	Value
Product Name	ITK/EMT Rabbit mAb
Remarks/Alias	EMT; LYK; LPFS1; PSCTK2; ITK/EMT
Species	Human
GeneID (Human)	3702
Immunogen	Recombinant protein of human ITK/EMT
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q08881
Protein Weight	72kDa
Shipping	Ice bag