


Anti-Itln1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wwa0

Item Number: CM0028742



Introduction

Rabbit polyclonal antibody against human ITLN1 (Intelectin-1/Omentin) validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat; ideal for innate immunity and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	ITLN1 Rabbit pAb
Remarks/Alias	HL1; LFR; HL-1; INTL; ITLN; hIntL; omentin; ITLN1
Species	Human
GeneID (human)	55600
GeneID	55600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-313 of human ITLN1 (NP_060095.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WWA0
Protein Weight	36kDa
Shipping	Ice bag