


Anti-Itpa Polyclonal Antibody For Wb, Elisa - Q9By32

Item Number: CM0016040



Introduction

Rabbit polyclonal antibody against human ITPA, a nucleotide metabolism enzyme. Validated for WB and ELISA; cross-reacts with human, mouse, rat. Suitable for studying non-canonical purine exclusion.

[Learn More](#)

Parameter	Specification
Product Name	ITPA Rabbit pAb
Remarks/Alias	DEE35; My049; ITPase; NTPase; C20orf37; dj794i6.3; HLC14-06-P; ITPA
Species	Human
Gene ID (Human)	3704
Gene ID	3704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human ITPA (NP_258412.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY32
Protein Weight	21kDa
Shipping	Ice bag