


Anti-Itpka Rabbit Polyclonal Antibody For Wb, Elisa - P23677

Item Number: CM0021942



Introduction

Rabbit polyclonal antibody targeting human ITPKA (inositol-trisphosphate 3-kinase A), suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying InsP3 metabolism and calcium signaling.

[Learn More](#)

Parameter	Value
Product Name	ITPKA Rabbit pAb
Remarks/Alias	IP3KA, IP3-3KA, ITPKA
Species	Human
Gene ID (Human)	3706
Gene ID	3706
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 141-461 of human ITPKA (NP_002211.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P23677
Protein Weight	51kDa
Shipping	Ice bag