

Anti-Itpkb Rabbit Polyclonal Antibody For Wb And Elisa - P27987

Item Number: CM0025943

Introduction

Rabbit polyclonal antibody against human ITPKB (IP3K-B), validated for WB and ELISA, cross-reactive with mouse and rat. Ideal for studying calcium homeostasis and inositol phosphate signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	ITPKB Rabbit pAb
Remarks/Alias	IP3K; IP3KB; PIG37; IP3K-B; IP3-3KB; ITPKB
Species	Human
GeneID (Human)	3707
GeneID	3707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ITPKB (NP_002212.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27987
Protein Weight	102kDa
Shipping	Ice bag