

Anti-Itpr2 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q14571

Item Number: CM0023514

Introduction

Rabbit polyclonal antibody targeting human ITPR2 (IP3R2), an IP3-gated calcium channel. Validated for IHC-P, IF/ICC, ELISA; cross-reacts with mouse, rat. Suitable for calcium signaling and muscle/heart research.

[Learn More](#)

Parameter	Specification
Product Name	ITPR2 Rabbit pAb
Remarks/Alias	ANHD; IP3R2; CFAP48; INSP3R2; ITPR2
Species	Human
Gene ID	3709
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14571
Protein Weight	308kDa
Shipping	Ice bag