


Anti-Itpr3 Rabbit Monoclonal Antibody For Wb, Elisa - Q14573

Item Number: CM0004035



Introduction

Rabbit monoclonal antibody targeting human ITPR3 (IP3R3) for Western blotting and ELISA. Ideal for studying inositol 1,4,5-trisphosphate-mediated calcium release and intracellular signaling.

[Learn More](#)

Parameter	Value
Product Name	ITPR3 Rabbit mAb
Remarks/Alias	IP3R; CMT1J; IP3R3; ITPR3
Species	Human
GeneID (Human)	3710
GeneID	3710
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1630-1870 of human ITPR3 (NP_002215.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14573
Protein Weight	304kDa
Shipping	Ice bag