


Anti-Itpr3 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q14573

Item Number: CM0018363



Introduction

High-quality rabbit polyclonal antibody targeting ITPR3 (IP3R-3), a calcium channel. Suitable for WB, IF/ICC, IP, ELISA on human and mouse samples. Ideal for cell signaling and calcium homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	ITPR3 Rabbit pAb
Remarks/Alias	IP3R; CMT1J; IP3R3; ITPR3
Species	Human
GeneID (Human)	3710
GeneID	3710
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1630-1870 of human ITPR3 (NP_002215.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14573
Protein Weight	304kDa
Shipping	Ice bag