

Anti-Itprrip Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iwb1

Item Number: CM0017157

Introduction

Rabbit polyclonal antibody against human ITPRIP (DANGER), specific for WB, IF/ICC, ELISA. Ideal for calcium signaling studies involving inositol 1,4,5-trisphosphate receptor interaction and membrane biology.

[Learn More](#)

Parameter	Value
Product Name	ITPRIP Rabbit pAb
Remarks/Alias	DANGER; KIAA1754; bA127L20; bA127L20.2; ITPRIP
Species	Human
GeneID (Human)	85450
GeneID	85450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-290 of human ITPRIP (NP_203755.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8IWB1
Protein Weight	62kDa
Shipping	Ice bag