


Anti-Ivd Polyclonal Antibody For Wb, Ihc-P, Elisa - P26440

Item Number: CM0023115



Introduction

Rabbit polyclonal antibody targeting human IVD (isovaleryl-CoA dehydrogenase). Cross-reacts with mouse and rat. Suitable for Western blot, IHC-P, and ELISA. Useful for studies of leucine catabolism and mitochondrial fatty acid oxidation.

[Learn More](#)

Parameter	Value
Product Name	IVD Rabbit pAb
Remarks/Alias	IVDH; ACAD2; IVD
Species	Human
Gene ID (Human)	3712
Gene ID	3712
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-200 of human IVD (NP_002216.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26440
Protein Weight	47kDa
Shipping	Ice bag